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## Caries arrest, remineralization and desensitization effect

### Silver Diamine Fluoride

FAgamin® contains diamine silver fluoride (Ag (NH<sub>3</sub>)<sub>2</sub>F). Other similar products used in caries prevention and dental hypersensitivity have several disadvantages in opposition to FAgamin®. When applying silver nitrate (AgNO<sub>3</sub>) or sodium fluoride (NaF), the dental enamel hydroxiapatite dissociates resulting in soluble calcium and phosphate ions. The action on teeth of diamine silver fluoride is based on the chemical reaction with hydroxiapatite.

Calcium fluoride (CaF<sub>2</sub>) and silver phosphate (Ag<sub>3</sub> PO<sub>4</sub>) form deposit over the tooth surface, reinforcing and protecting it. This explains the remineralizable effect of FAgamin®.

Scientific studies reveal that silver ion (Ag<sup>+</sup>) inhibits bacteria growth and has a strong action against enzymes. Fluoride ion (F<sup>-</sup>) reduces hydroxiapatite solubility and inhibits caries proliferation.

In case of applying FAgamin® over dental caries, the result will be a dark deposit due to the reaction of silver ions with the organic material. This is why it can also play the role of a caries detector.

FAgamin® possesses the advantage of having a wide range of clinical applications attributable to its chemical composition. Some of these are: old teeth enamel and caries nearby dentin remineralization, caries detection and arrest, incipient and old caries inhibition, dental hypersensitivity desensitization. It also has an antibacterial and antiseptic action.

Furthermore, FAgamin® offers a practical, easy and quick topical method for dental caries treatment on temporary teeth.

#### Applications:

It can be applied on cavities prepared for prosthesis, crowns, bridges, incipient or rampant dental caries and in fissures and cavities (after widening) if dentist considers so. It can also be applied in infected root canals, diluted 1:10 in distilled water.

#### Instructions for use:

Wash and dry the tooth before beginning the treatment.

Damp a cotton swab or soft brush with two or three drops of FAgamin®, wait a few minutes and rub thoroughly over the tooth for one or two minutes. When finished, wash to remove the solution residue.

To avoid the contact of FAgamin® with the mucous membrane, isolate it with a cotton pellet or rubber dump.

In cases where the carious tissue is near gingival tissue, cover it with Vaseline or cacao butter.

If, accidentally, the contact occurs or the patient shows pain, immediately wash water with salt or oxygenated water (10 volumes).

**Warnings:**

- As the reaction of silver could darken the dentin, it is not convenient to use FAgamin® on the frontal permanent teeth.
- In deep dental caries cavities, as a filtration to the pulp is likely to occur (due to the tooth quality for instance), it is advisable to avoid its use or dilute 1:10 with distilled water.

**Precautions:**

- Use gloves.
- Attach the bottle cap immediately after using.
- Store in a place protected from light.
- Not to be stored at temperatures exceeding 20°C.
- If contact with skin occurs, wash with ammonium solution or oxygenated water.

**Contraindications:**

This product has no contraindications in the doses to be applied.

**Composition:**

8% Diamine Silver Fluoride solution  $\text{Ag}(\text{NH}_3)_2\text{F}$ .

**Content:**

- FAgamin® contains 5 ml of 38% Diamine Silver Fluoride solution. Controlled pH.

**Storage:**

Store in a dark, cool (between 4 and 30 °C) and dry place. Keep the container tightly closed.

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